



Point-Of-Use Performance Chart



Hydro Picosystems are designed to remove all contaminants, all of the time. Picosystems can use source water from the tap, or from larger ultrapure water systems, depending on your application's specific needs. Our Picosystems are designed for flexibility, so that you can choose the system that meets your specific requirements. When your applications change, you can easily add components to upgrade your system. Your Picosystem can be conveniently stored on a shelf, under a counter, or wherever your lab space permits.

	PICOSYSTEM	PICOSYSTEM PLUS	PICOSYSTEM UV PLUS	PICOTECH 2	PICOTECH 2 PLUS	PICOTECH 2 UV PLUS	PICOPURE 2	PICOPURE 2 PLUS	PICOPURE 2 UV PLUS
18 megohm-cm resistivity	●	●	●	●	●	●	●	●	●
< 5 ppb TOC			●			●			●
< 20 ppb TOC		●			●			●	
< 35 ppb TOC	●			●			●		
< 0.1 CFU/mL bacteria				●	●	●	●	●	●
< 0.1 CFU/mL bacteria with capsule post filter	●	●	●						
< 5 CFU/mL bacteria			●						
< 10 CFU/mL bacteria	●	●							
endotoxin & pyrogen free, < 0.25 EU/mL							●	●	●
RNase free up to 20 hours									●